

Technical Specification

PRODUCT : IODINE 0.1 N (0.05 mol/L) RPE - For analysis
CODE : 456030000

TEST	U.M.	SPECIFICATION
Description	-	Brown red liquid
Assay (potentiometry)	N	0.0998 - 0.1002

This volumetric solution is checked by sodium thiosulfate solution, tested with potassium dichromate volumetric standard traceable to the SRM 136f of NIST (National Institute of Standards and Technology)

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in